

FIRM LOCK PLATES

Firm-lock T-buttress plate 4.5/5.0
 Firm-lock placas de sostén 4.5/5.0 en T
 Firm-lock plaques de soutien 4.5/5.0 en T



Stainless Steel	Size	Length mm	Titanium
Item Code			Item Code
626.004	4 Holes	82	626T.004
626.005	5 Holes	98	626T.005
626.006	6 Holes	114	626T.006
626.007	7 Holes	130	626T.007
626.008	8 Holes	146	626T.008
626.009	9 Holes	162	626T.009
626.010	10 Holes	178	626T.010
626.011	11 Holes	194	626T.011
626.012	12 Holes	210	626T.012
626.013	13 Holes	226	626T.013
626.014	14 Holes	242	626T.014
626.016	16 Holes	274	626T.016

*Without LC Undercuts. All holes in shaft counted.
 Sizes not listed above can also be manufactured on request.*