


FIRM LOCK PLATES

Firm-lock distal dorsal distal radius L-plate 2.4 - right angled

Firm-lock placas dorsales, 2.4 en L para radio distal, acodada hacia la derecha

Firm-lock plaques dorsales, 2.4 en L pour radius distal, coudée à droite



Stainless Steel		Head	Shaft	Length	Titanium	
Item Code					Item Code	
left angled	right angled				left angled	right angled
671.003L	671.003R	2 Holes	3Holes	40	671T.003L	671T.003R
671.004L	671.004R	2 Holes	4Holes	49	671T.004L	671T.004R
671.005L	671.005R	2 Holes	5Holes	58	671T.005L	671T.005R

*With LC undercuts. All holes in shaft counted.
 Sizes not listed above can also be manufactured on request.*