

FIRM LOCK PLATES

Firm-lock superior anterior clavicle plate 2.7/3.5 with lateral extension



Stainless Steel		Size	Length mm	Titanium	
Item Code				Item Code	
left	right			left	right
90586.003	90585.003	3 Holes	69	90586T.003	90585T.003
90586.004	90585.004	4 Holes	81	90586T.004	90585T.004
90586.005	90585.005	5 Holes	94	90586T.005	90585T.005
90586.006	90585.006	6 Holes	108	90586T.006	90585T.006
90586.007	90585.007	7 Holes	123	90586T.007	90585T.007
90586.008	90585.008	8 Holes	135	90586T.008	90585T.008

Sizes not listed above can also be manufactured on request.