

# ETO STERILIZER



## Specifications :

- Sterility assurance through automatic temperature control ( $\pm 3^{\circ}\text{C}$ )
- Accurate control of temperature, moisture vacuum valves and all critical parameters.
- PLC based automation with 6 pre-programmed cycles & additional user defined cycles.
- Automatic aeration cycle at the end of process for ETO exhaust.
- Operates on universally available 100% ETO Cartridge.
- Live display of phase status, current value, set values, temperature, pressure & humidity.
- Current & previous reports can be printed through an external printer / internal printer.

## Safety Features :

- ETO cartridge punctured only in negative pressure & humidity conditions.
- Chamber fabricated with heavy duty SS 304 / SS 316.
- Glass Wool insulation.

Volume Ltrs	25	56	85	110	190	256	450
Internal Dimension (inch)	12x12x12	12x12x24	12x12x36	12x12x48	18x18x36	18x18x48	24x24x48